



Release Note	
Product Name: Pixyl.Neuro	Release Date: 31-Aug-2020

Release Note for Pixyl.Neuro

Version Number
Version 1.6.2
Introduction
Pixyl.Neuro is a web-based software to help diagnosis in neurology and neuroradiology based on brain MRI data. Pixyl.Neuro performs biomarkers extraction for multiple sclerosis, pathologies that affect the white matter in general (Leukoencephalopathy for example), and neurodegenerative diseases (Alzheimer's, Parkinson, etc.). It provides fully automatic brain lesion segmentation, cortical and subcortical structures segmentation, comparison with a normal population, and longitudinal analysis.
Back-end release notes
1.6.2 (2020-08-31) <ul style="list-style-type: none">Improved the calculation of the new lesions found in the longitudinal MS/WMH analysisFixed the reproducibility for analysis running on GPUsKeep the DICOM tag "SeriesDate" in segmented FLAIR
Front-end release notes
1.6.2 (2020-08-31) <ul style="list-style-type: none">Management of the new Compute Engine resource creation rate limits enforced by the Cloud provider
Release Authorised by
Senan Doyle, CEO