



Release Note	
Product Name: Pixyl.Neuro	Release Date: 09-Nov-2020

Release Note for Pixyl.Neuro

Version Number
Version 1.6.3
Introduction
Pixyl.Neuro is a web-based software to help diagnosis in neurology and neuroradiology based on brain MRI data. Pixyl.Neuro performs biomarkers extraction for multiple sclerosis, pathologies that affect the white matter in general (Leukoencephalopathy for example), and neurodegenerative diseases (Alzheimer's, Parkinson, etc.). It provides fully automatic brain lesion segmentation, cortical and subcortical structures segmentation, comparison with a normal population, and longitudinal analysis.
Back-end release notes
<ul style="list-style-type: none">•
Front-end release notes
1.6.3 (2020-11-09) <ul style="list-style-type: none">• Fixed a potential risk of data mix between two analyses launched in parallel
Release Authorised by
Senan Doyle, CEO